



PRODUCT PROFILE

- PRODUCT:** Purified Mouse anti-FeLV (Clone #8D11)
- CATALOG NUMBER:** MFLV-20M-8D11
- LOT NUMBER:**
- QUANTITY:**
- PHYSICAL STATE:** Liquid
- BUFFER:** PBS, pH 7.2
- PRESERVATIVE:** 0.09% Sodium Azide
- STORAGE:** 2 - 8° C for less than 30 days. Aliquot and store between -20C to -80C for long term storage. Avoid multiple freeze-thaw cycles.
- EXPIRATION:** 12 months from shipping date.
- PROCEDURE:** This mouse IgG1 kappa is purified using salt fractionation procedures and Protein A chromatography. It is produced from tissue culture media obtained at our facility in The United States. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.
- SPECIFICATIONS:** This antibody reacts with feline leukemia virus (FeLV) p27 as determined by ELISA, LFD and western blot. Suggest using this as ELISA detection antibody when capturing with MFLV-20M-CM1. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR CAT OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA