



PRODUCT PROFILE

PRODUCT: Purified Mouse anti-Giardia Lamblia (Clone #4G5)

CATALOG NUMBER: MGRD-55A-4G5

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: PBS, pH 7.2

PRESERVATIVE: 0.09% Sodium Azide

STORAGE: 2 - 8° C

EXPIRATION: 12 months from shipping date.

PROCEDURE: This mouse IgG1 kappa is purified using salt fractionation procedures and Protein A chromatography. It is produced from tissue culture media obtained at our facility in The United States. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.

SPECIFICATIONS: This antibody reacts with cysts from giardia lamblia as determined by blotting, ELISA, IF and LFD techniques. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA