



PRODUCT PROFILE

PRODUCT: Purified Mouse anti-Plasmodium falciparum HRP11

CATALOG NUMBER: MPFM-55A

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: PBS, pH 7.2

PRESERVATIVE: 0.09% Sodium Azide

STORAGE: 2 - 8° C

EXPIRATION: 12 months from shipping date.

PROCEDURE: This mouse IgM kappa is purified using salt fractionation procedures and column chromatography. Antibody concentration was determined using an absorbance at 280 nm: 1.2 equals 1.0 mg of IgM.

SPECIFICATIONS: This antibody reacts with native histidine rich protein 2 (HRP11) from Plasmodium falciparum as determined by blotting assays and ELISA techniques. It is shown to pair with ICL's MPFG-55A in sandwich ELISA and Lateral Flow. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA