

Immunology Consultants Laboratory, Inc.

7150 SW Sandburg Street
Portland, OR 97223
T 503-747-2454
F 503-747-2544

www.icllab.com



ANTIBODIES ARE OUR BUSINESS

PRODUCT PROFILE

PRODUCT: Purified Mouse anti-Human TSH beta (Clone #3B9)

CATALOG NUMBER: MTSH-80A-3B9

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: PBS, pH 7.2.

PRESERVATIVE: 0.09% sodium azide

STORAGE: Antibody should be stored between 2 – 8 °C for up to 14 days. For long-term storage, aliquot and store between -20°C and -80°C. Multiple freeze-thaw cycles are not recommended.

EXPIRATION: 12 months from shipping date.

PROCEDURE: This mouse IgG1 kappa is purified using salt fractionation procedures and Protein A chromatography. It is produced from tissue culture media obtained at our facility in The United States. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.

SPECIFICATIONS: This antibody reacts with native human thyroid stimulating hormone beta subunit (TSHb) as determined by ELISA and western blot. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA

MTSH-80A-3B9_Rev030119