



PRODUCT PROFILE

PRODUCT: Purified Mouse anti-Zika Virus NS1 (Clone #2B9)

CATALOG NUMBER: MZVN-55A-2B9

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: PBS, pH 7.2

PRESERVATIVE: 0.09% sodium azide

STORAGE: 2 – 8 C

EXPIRATION: 12 months from shipping date.

PROCEDURE: This mouse IgG1 kappa is purified using salt fractionation procedures and Protein A chromatography. It is produced from tissue culture media obtained at our facility in The United States. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.

SPECIFICATIONS: This antibody reacts with recombinant Zika Virus NS1 of as determined by ELISA and blotting techniques. Optimal working dilutions should be determined for the particular assay conditions.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA