

Immunology Consultants Laboratory, Inc.

7150 SW Sandburg Street
Portland, OR 97223
T 503-747-2454
F 503-747-2544

www.icllab.com



ANTIBODIES ARE OUR BUSINESS

PRODUCT PROFILE

PRODUCT: Purified Chicken IgY Isotype Control

CATALOG NUMBER: RA-30Y

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: Phosphate buffered saline, pH 7.2

PRESERVATIVE: 0.09% sodium azide

STORAGE: 2 -8°C

PROCEDURE: The IgY was purified by a combination of extraction and salt fractionation.

SPECIFICATIONS: Antibody concentration was determined by extinction coefficient : O.D. 1.32 at 280nm equals 1.0 mg of IgY (IgG). Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA

RA30Y_Rev061919