

Immunology Consultants Laboratory, Inc.

7150 SW Sandburg Street
Portland, OR 97223
T 503-747-2454
F 503-747-2544

www.icllab.com



ANTIBODIES ARE OUR BUSINESS

PRODUCT PROFILE

PRODUCT: Dog Heartworm (*D. immitis*) Antigen

CATALOG NUMBER: RA-40HRT

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: PBS, pH 7.2

PRESERVATIVE: 0.09% sodium azide

STORAGE: 2 - 8° C for less than 30 days. Aliquot and store between -20C to -80C for long term storage. Avoid multiple freeze-thaw cycles.

EXPIRATION: 12 months from shipping date.

PROCEDURE: Washed, adult female worms were disrupted, homogenized and protein was extracted. Total protein measurement was determined by Bio-Rad Protein Assay using BGG as a standard.

SPECIFICATIONS: This product shows multiple bands as determined by SDS-PAGE gel. It is suitable for use as a control in ELISA, blotting or other immunoassays. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA