



PRODUCT PROFILE

PRODUCT: Affinity Purified Rabbit anti-BSA

CATALOG NUMBER: RAL-10A

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: PBS, pH 7.2

PRESERVATIVE: 0.09% sodium azide

STORAGE: 2 - 8° C

EXPIRATION: 24 months from shipping date.

PROCEDURE: Rabbits were immunized with highly purified bovine serum albumin (BSA) and the resulting antiserum was collected. Antibody was immunoaffinity purified off an albumin containing immunosorbent. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.

SPECIFICATIONS: This antibody reacts with BSA as determined by IEP and ELISA techniques. It is suitable for blotting, ELISA and lateral flow. This antibody may cross react with albumin from other species. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR BOVINE OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA