



PRODUCT PROFILE

PRODUCT: HRP Conjugated Rabbit anti-CBP Tag

CATALOG NUMBER: RCBP-45P

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: HEPES – saline with stabilizer

PRESERVATIVE: 0.002% thimerosal

STORAGE: 2 - 8° C

EXPIRATION: 12 months from shipping date.

PROCEDURE: Rabbits were immunized with highly purified MKRRWKKNFIAVSAANRFKKISSSGAL, Calmodulin Binding Protein (CBP) tag, conjugated with KLH and the resulting antiserum was collected. Antibody was immunoaffinity purified off a peptide containing immunosorbent and conjugated with horseradish peroxidase (HRP). Antibody concentration was determined prior to conjugation using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.

SPECIFICATIONS: This antibody reacts with MKRRWKKNFIAVSAANRFKKISSSGAL as determined by IEP and ELISA techniques. It is suitable for blotting, ELISA and IHC applications. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA