



PRODUCT PROFILE

PRODUCT: Rabbit anti-Dog CRP

CATALOG NUMBER: RCRP-40

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: None added.

PRESERVATIVE: 0.09% sodium azide

STORAGE: 2 - 8° C

EXPIRATION: 24 months from shipping date.

PROCEDURE: Rabbits were immunized with highly purified dog c-reactive protein (CRP) and the resulting antiserum was collected. Total protein measurement was determined by refractometer. Specific antibody content is determined immunoaffinity purification off an antigen containing immunosorbent.

SPECIFICATIONS: This antibody reacts with dog CRP as determined by IEP and ELISA techniques. It is suitable for blotting, ELISA and various other immunoassay platforms. This antibody may cross react with CRP from other species. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA