



PRODUCT PROFILE

PRODUCT: HRP Conjugated Rabbit anti S-tag

CATALOG NUMBER: RES-45P-Z

LOT NUMBER:

QUANTITY:

BUFFER: HEPES – saline with stabilizer

PRESERVATIVE: 0.002% thimerosal

STORAGE: 2 - 8° C

EXPIRATION: 12 months from shipping date.

PROCEDURE: Rabbits were immunized with highly purified S-tag peptide conjugated with KLH and the resulting antiserum was collected. Antibody was solid phase adsorbed, immunoaffinity purified off a peptide containing immunosorbent and conjugated with horseradish peroxidase (HRP). Antibody concentration was determined prior to conjugation using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.

SPECIFICATIONS: This antibody reacts with S-tag as determined by IEP and ELISA techniques. It is suitable for blotting, ELISA and IHC applications. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA