

# Immunology Consultants Laboratory, Inc.

7150 SW Sandburg Street  
Portland, OR 97223  
T 503-747-2454  
F 503-747-2544

[www.icllab.com](http://www.icllab.com)



ANTIBODIES ARE OUR BUSINESS

## PRODUCT PROFILE

**PRODUCT:** HRP Conjugated Rabbit anti S-tag

**CATALOG NUMBER:** RES-45P

**LOT NUMBER:**

**QUANTITY:**

**BUFFER:** HEPES – saline with stabilizer

**PRESERVATIVE:** 0.002% thimerosal

**STORAGE:** 2 - 8° C

**EXPIRATION:** 12 months from shipping date.

**PROCEDURE:** Rabbits were immunized with highly purified S-tag peptide conjugated with KLH and the resulting antiserum was collected. Antibody was solid phase adsorbed, immunoaffinity purified off a peptide containing immunosorbent and conjugated with horseradish peroxidase (HRP). Antibody concentration was determined prior to conjugation using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.

**SPECIFICATIONS:** This antibody reacts with S-tag as determined by IEP and ELISA techniques. It is suitable for blotting, ELISA and IHC applications. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.  
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

***PRODUCT OF THE UNITED STATES OF AMERICA***