



PRODUCT PROFILE

PRODUCT:	Affinity Purified Rabbit anti-beta Galactosidase
CATALOG NUMBER:	RGAL-45A-Z
LOT NUMBER:	
QUANTITY:	1.0 mL @ 1.0 mg/mL
PHYSICAL STATE:	Liquid
BUFFER:	PBS, pH 7.2
PRESERVATIVE:	0.09% sodium azide
STORAGE:	2 - 8° C
EXPIRATION:	24 months from shipping date.
PROCEDURE:	Rabbits were immunized with highly purified beta-galactosidase from <i>Escherichia coli</i> . and the resulting antiserum was collected. Antibody was immunoaffinity purified off a peptide containing immunosorbent. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.
SPECIFICATIONS:	This antibody reacts with beta-galactosidase as determined by IEP and ELISA techniques. It is suitable for blotting, ELISA and IHC applications. This antibody may cross react with beta-galactosidase from other species. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA