



PRODUCT PROFILE

PRODUCT: Affinity Purified Rabbit anti-Bovine IgG h+l

CATALOG NUMBER: RGHL-10A

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: PBS, pH 7.2

PRESERVATIVE: 0.09% sodium azide

STORAGE: 2 - 8° C

EXPIRATION: 24 months from shipping date.

PROCEDURE: Rabbits were immunized with highly purified bovine IgG h+l and the resulting antiserum was collected. Specific antibody was immunoaffinity purified off a bovine IgG containing immunosorbent. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.

SPECIFICATIONS: This antibody will react with bovine IgG and with light chains common with other bovine immunoglobulins as determined by ELISA and IEP techniques. It is suitable for blotting, ELISA, IHC and Lateral Flow applications. This antibody may cross-react with IgG from other species. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR BOVINE OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA