



## **PRODUCT PROFILE**

**PRODUCT:** Affinity Purified Rabbit anti-Sheep IgG h+l

**CATALOG NUMBER:** RGHL-35A

**LOT NUMBER:**

**QUANTITY:**

**PHYSICAL STATE:** Liquid

**BUFFER:** PBS, pH 7.2

**PRESERVATIVE:** 0.09% sodium azide

**STORAGE:** 2 - 8° C

**EXPIRATION:** 24 months from shipping date.

**PROCEDURE:** Rabbits were immunized with highly purified sheep IgG h+l and the resulting antiserum was collected. Specific antibody was immunoaffinity purified off a sheep IgG containing immunosorbent. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.

**SPECIFICATIONS:** This antibody will react with sheep IgG and with light chains common with other sheep immunoglobulins as determined by ELISA and IEP techniques. It is suitable for blotting, ELISA, IHC and Lateral Flow applications. This antibody may cross-react with IgG from other species. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.  
FOR IN VITRO USE ONLY, NOT FOR SHEEP OR ANIMAL CONSUMPTION**

**PRODUCT OF THE UNITED STATES OF AMERICA**