

Immunology Consultants Laboratory, Inc.

7150 SW Sandburg Street
Portland, OR 97223

T 503-747-2454

F 503-747-2544

www.icllab.com



ANTIBODIES ARE OUR BUSINESS

PRODUCT PROFILE

PRODUCT: FITC Conjugated Rabbit anti-Goat IgG h+l

CATALOG NUMBER: RGHL-90F

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: PBS, pH 7.2 with 0.2% BSA.

PRESERVATIVE: 0.09% sodium azide

STORAGE: 2 - 8° C

EXPIRATION: 12 months from shipping date.

PROCEDURE: Rabbits were immunized with highly purified goat IgG h+l and the resulting antiserum was collected. Specific antibody was immunoaffinity purified off a goat IgG containing immunosorbent and conjugated with fluorescein isothiocyanate (FITC). Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG. FITC/Protein Ratio = 3.0 – 6.0.

SPECIFICATIONS: This antibody will react with goat IgG and with light chains common with other goat immunoglobulins as determined by ELISA and IEP techniques. It is suitable for blotting, ELISA and IF applications. This antibody may cross-react with IgG from other species. Optimal working dilutions should be determined experimentally by the investigator.

FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION

PRODUCT OF THE UNITED STATES OF AMERICA