



PRODUCT PROFILE

- PRODUCT:** Agarose Immobilized Rabbit anti-HA Tag
- CATALOG NUMBER:** RHGT-145D
- LOT NUMBER:**
- QUANTITY:**
- PHYSICAL STATE:** Liquid Suspension (~50% slurry)
- BUFFER:** PBS, pH 7.2
- PRESERVATIVE:** 0.09% sodium azide
- STORAGE:** 2 - 8° C
- EXPIRATION:** 12 months after shipment date.
- PROCEDURE:** Rabbits were immunized with highly purified YPYDVPDYA (influenza hemagglutinin-HA-epitope) and the resulting antiserum was collected. Antibody was immunoaffinity purified off a peptide containing immunosorbent then coupled to agarose beads. Antibody concentration was determined prior to agarose coupling using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.
- SPECIFICATIONS:** This antibody reacts with YPYDVPDYA as determined by IEP and ELISA techniques. It is suitable for immunoprecipitation. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA