



## **PRODUCT PROFILE**

**PRODUCT:** Affinity Purified Rabbit anti-C-MYC

**CATALOG NUMBER:** RMYC-45A-Z

**LOT NUMBER:**

**QUANTITY:**

**PHYSICAL STATE:** Liquid

**BUFFER:** PBS, pH 7.2

**PRESERVATIVE:** 0.09% sodium azide

**STORAGE:** 2 - 8° C

**EXPIRATION:** 24 months from shipping date.

**PROCEDURE:** Rabbits were immunized with highly purified EQKLISEEDL conjugated with KLH (C-MYC) and the resulting antiserum was collected. Antibody was solid phase adsorbed then immunoaffinity purified off an antigen containing immunosorbent. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.

**SPECIFICATIONS:** This antibody reacts with EQKLISEEDL as determined by IEP and ELISA techniques. It is suitable for blotting, ELISA and IHC applications. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.  
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

**PRODUCT OF THE UNITED STATES OF AMERICA**