



PRODUCT PROFILE

PRODUCT:	Affinity Purified Rabbit anti-Streptavidin
CATALOG NUMBER:	RSVD-65A
LOT NUMBER:	
QUANTITY:	
PHYSICAL STATE:	Liquid
BUFFER:	PBS, pH 7.2
PRESERVATIVE:	0.09% Sodium azide
STORAGE:	2 – 8 C
EXPIRATION:	24 months from shipping date.
PROCEDURE:	Rabbits were immunized with Streptavidin and the resulting antiserum was collected. Antibody was immunoaffinity purified off a Streptavidin containing immunosorbent. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.
SPECIFICATIONS:	This antibody reacts with Streptavidin as determined by IEP and ELISA techniques. It is suitable for blotting and ELISA applications. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA