



PRODUCT PROFILE

PRODUCT: HRP Conjugated Rabbit anti-T7 tag

CATALOG NUMBER: RT7-45P-Z

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: HEPES – saline with stabilizer

PRESERVATIVE: 0.002% thimerosal

STORAGE: 2 - 8° C

EXPIRATION: 12 months from shipping date.

PROCEDURE: Rabbits were immunized with highly purified T7 epitope tag conjugated to KLH. and the resulting antiserum was collected. Antibody was solid phase adsorbed, immunoaffinity purified off a peptide containing immunosorbent and conjugated with horseradish peroxidase (HRP). Antibody concentration was determined prior to conjugation using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.

SPECIFICATIONS: This antibody reacts with T7 as determined by IEP and ELISA techniques. It is suitable for blotting, ELISA and IHC applications. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA