



PRODUCT PROFILE

PRODUCT: FITC Conjugated Rabbit anti-V5 tag

CATALOG NUMBER: RV5-45F-Z

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: PBS, pH 7.2 with 0.2% BSA.

PRESERVATIVE: 0.09% sodium azide

STORAGE: 2 - 8° C

EXPIRATION: 12 months from shipping date.

PROCEDURE: Rabbits were immunized with highly purified V5 peptide conjugated with KLH and the resulting antiserum was collected. Specific antibody was immunoaffinity purified off a peptide containing immunosorbent and conjugated with fluorescein isothiocyanate (FITC). Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG. FITC/Protein Ratio = 3.0 – 6.0.

SPECIFICATIONS: This antibody will react with V5 by ELISA and IEP techniques. It is suitable for blotting, ELISA and IF applications. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA