

Immunology Consultants Laboratory, Inc.

7150 SW Sandburg Street
Portland, OR 97223
T 503-747-2454
F 503-747-2544

www.icllab.com



ANTIBODIES ARE OUR BUSINESS

PRODUCT PROFILE

PRODUCT: Rabbit anti-Bovine Whole Serum

CATALOG NUMBER: RWS-10

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: N/A

PRESERVATIVE: 0.09% sodium azide

STORAGE: 2 - 8° C

EXPIRATION: 24 months from shipping date.

PROCEDURE: Rabbits were immunized with normal bovine serum and the resulting antiserum was collected. Antibody concentration was determined using a refractometer.

SPECIFICATIONS: This antibody reacts with multiple proteins in bovine serum when tested by immunoelectrophoresis (IEP). Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR BOVINE OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA

RWS-10_Rev031319