



PRODUCT PROFILE

PRODUCT: Rabbit anti-Cat Whole Serum

CATALOG NUMBER: RWS-20

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: N/A

PRESERVATIVE: 0.09% sodium azide

STORAGE: 2 - 8° C

EXPIRATION: 24 months from shipping date.

PROCEDURE: Rabbits were immunized with normal cat serum and the resulting antiserum was collected. Antibody concentration was determined using a refractometer.

SPECIFICATIONS: This antibody reacts with multiple proteins in cat serum when tested by immunoelectrophoresis (IEP). Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA