



PRODUCT PROFILE

PRODUCT:	Rabbit anti-Pig Whole Serum
CATALOG NUMBER:	RWS-5
LOT NUMBER:	
QUANTITY:	
PHYSICAL STATE:	Liquid
BUFFER:	N/A
PRESERVATIVE:	0.09% sodium azide
STORAGE:	2 - 8° C
EXPIRATION:	24 months from shipping date.
PROCEDURE:	Rabbits were immunized with normal pig serum and the resulting antiserum was collected. Antibody concentration was determined using a refractometer.
SPECIFICATIONS:	This antibody reacts with multiple proteins in pig serum when tested by immunoelectrophoresis (IEP). Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA