

# Immunology Consultants Laboratory, Inc.

7150 SW Sandburg Street  
Portland, OR 97223  
T 503-747-2454  
F 503-747-2544

[www.icllab.com](http://www.icllab.com)



ANTIBODIES ARE OUR BUSINESS

## PRODUCT PROFILE

**PRODUCT:** Affinity Purified Sheep anti-Dog CRP

**CATALOG NUMBER:** SCRP-40A

**LOT NUMBER:**

**QUANTITY:**

**PHYSICAL STATE:** Liquid

**BUFFER:** PBS, pH 7.2

**PRESERVATIVE:** 0.09% sodium azide

**STORAGE:** 2 - 8° C

**EXPIRATION:** 24 months from shipping date.

**PROCEDURE:** Sheep were immunized with highly purified dog c-reactive protein (CRP) and the resulting antiserum was collected. Antibody was immunoaffinity purified off a dog CRP containing immunosorbent. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.

**SPECIFICATIONS:** This antibody reacts with dog CRP as determined by IEP and ELISA techniques. It is suitable for blotting, ELISA and IHC applications. This antibody may cross react with CRP from other species. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.  
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

**PRODUCT OF THE UNITED STATES OF AMERICA**

SCRP-40A\_Rev031219