

## **PRODUCT PROFILE**

**PRODUCT:** Goat anti-Rabbit Whole Serum

**CATALOG NUMBER:** GWS-15

**LOT NUMBER:**

**QUANTITY:**

**PHYSICAL STATE:** Liquid

**BUFFER:** N/A

**PRESERVATIVE:** 0.09% sodium azide

**STORAGE:** 2 - 8° C

**EXPIRATION:** 24 months from shipping date.

**PROCEDURE:** Goats, born and raised in the USA, were immunized with normal rabbit serum. The resulting antiserum was collected, delipidated, and 0.45 µm filtered. Total protein is determined by using a refractometer.

**SPECIFICATIONS:** This antibody reacts with multiple proteins in rabbit serum when tested by immunoelectrophoresis (IEP). Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.  
FOR IN VITRO USE ONLY, NOT FOR RABBIT OR ANIMAL CONSUMPTION**

***PRODUCT OF THE UNITED STATES OF AMERICA***